

## **STERILIZATION: Tungsten Eye Shields and Superficial Silver Plated Lead Eye Shields**

Tungsten Eye Shields and Superficial Silver Plated Lead Eye Shields have been validated for use in the following STERRAD® Systems, manufactured by Advanced Sterilization Products, a Johnson & Johnson Company:

- STERRAD® 50 System
- STERRAD® 200 System
- STERRAD® 100S System
- STERRAD® NX™ System
- STERRAD® 100NX™ System

Below is a complete list of Medical Devices validated for use in STERRAD® Systems: STERRAD 50, 100S, 200, NX, 100NX:

Item Number	Description
936-583 to 627	Tungsten Eye Shields
934-012 to 018	Silver Plated Lead Eye Shields

Testing was performed in accordance with AAMI TIR No.12-2004 guidelines, “Designing, Testing and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers.”

Please refer to the STERRAD System User’s Guide for general reprocessing instructions, including proper cleaning, and drying, and packaging information prior to reprocessing any medical device in a STERRAD System. Specific sterilization questions should be made directly to Advanced Sterilization Products by calling 1-888-STERRAD.

All medical devices should be inspected prior to use to ensure patient safety. Do not use any device if there is a defect or damage noted during the inspection.