### J DENTISTRY





This product is available through:

**JRT** Associates 800-221-0111



25

## **Pro-Dent Set**

The Pro-Dent set is a **universal set of phantoms** for carrying out constancy and acceptance tests of **conventional and digital** dental X-ray units (intra-oral, panoramic and cephalometric).

This is **not an all-in-one device** where results of tests blur each other out. This is the **only solution on the market** that makes it possible to **measure the X-ray beam collimation** with a dental film or a digital detector.

The **Pro-Dent** set consists of:

## Pro-Dent $\alpha$ phantom

# Technical data (can be modified to customer specifications):

- 3 step wedge (first step made of copper foil 0.3 mm thick, the next ones are made of polytetrafluoroethylene 8 and 16 mm thick)
- cover made of PMMA

### Pro-Dent $\beta$ phantom

# Technical data (can be modified to customer specifications):

- cone for perpendicular X-Ray beam control in the range of  $0^\circ \div 1.5^\circ$
- pattern for beam radius estimation
- cover made of PMMA

universal set of phantoms for analogue and digital units

### Pro-Dent γ phantom

# Technical data (can be modified to customer specifications):

- pattern for line pair resolution evaluation (from 4 to 8 LP/mm)
- optional second pattern for line pair resolution evaluation (from 1.6 to 3 LP/mm)
  – OPG units
- four holes in 0.5 mm Al foil for low contrast resolution tests
- additional 6 mm aluminium filter
- cover made of PMMA

## } | DENTISTRY



This product is available through: **JRT** Associates

800-221-0111

#### With the Pro-Dent set you can do the following tests:

- dose reproducibility
- development process stability
- perpendicular X-ray beam (range  $0^{\circ} \div 1.5^{\circ}$ )
- limitation and alignment of the X-Ray beam (including beam radius measurement)
- spatial / line pair resolution (perpendicular, parallel and rotated 45° to anode-cathode line)
- low contrast resolution



26

QC



### **Product features:**

- complies with:
  - IEC 61223-2-7 and IEC 61223-3-4
  - DIN 6868-5 and DIN V 6868-151
  - ÖNORM S 5240-5 and ÖNORM S 5240-11
- CE certified
- the Manual provides detailed guidelines for carrying out each test,
- results assessment and registration

#### Accesories:

- 0.8 mm copper filter patient's head equivalent
- positioning stand with place for an analogue film or a digital detector
- band for firm, perpendicular attachment of the phantom to the X-Ray unit's beam applicator

nt Cu 0.8

E

- rings for centering phantoms on the X-Ray unit's beam applicator
- CD with documentation
- elegant and convenient aluminium box for storing phantoms



