





This product is available through:

JRT Associates 800-221-0111



Pro-RF Stand

Spacer plate with adjustable legs for positioning Pro-Alpha and Pro-Digi test phantoms near the X-Ray tube on units without insertion slots (e.g. over-table / under-table tubes, C-arms etc.).

Technical data (can be modified to customer specifications):

- table top made of PMMA
- legs made of aluminium
- · dedicated space for the Pro-RF Al 25 absorption filter

Product features:

- Complies with:IEC 61223-3-1
- CE certified
- the Manual provides detailed guidelines for carrying out each test, results assessment and registration























